

Y5 Series

Sleek modern designs and state of the art image quality are trademarks of the Y-Series cameras. The Y5 camera reaches a level of image quality normally reserved for still image photography of the highest caliber



The Y5 offers a class-leading combination of high resolution, high image quality and high speed. Starting with **2336 x 1728 pixels @ 730 fps (1455 fps Plus mode)** the Y5 produces the levels of resolution commercial applications demand combined with outstanding image quality rarely found in high speed video cameras.

High image quality using brand new innovations that produce ultra sharp images with very fine gradations even when shooting under less than ideal lighting conditions are part of the Y5 feature set. When the application demands faster frame rates, a unique region of interest or longer record times the Y5 has the capability. Typical of the Y-series cameras the Y5 offers fast Ethernet and USB2 connections with industry standard “save-as” digital image formats, LabView, MatLab drivers and fully integrated control software for MAC or PC.

Wireless Remote Control Option and more

- Full 2360 x 1728 @ 730 fps (1455 fps Plus mode)
- 30 to 125,000 user selectable pictures per second
- User selectable ROI
- 8GB or 16GB memory options
- I-Pod & I-Phone remote control option
- MAC and PC compatible software
- Gigabit Ethernet backward compatible with 100Base-T Ethernet
- Non-Volatile Image backup option
- LabView and MatLab Plug-ins.
- Sensitivity >2000 ASA mono, >800 ASA Color
- Size 90.5x90.5x220mm; Weight 1.9Kg

Contact: Digital West Imaging
450 Mountain View Road
San Diego, CA 92021

Sales@DigitalWestImaging.com
Telephone: 619-593-1900
FAX 619-593-1901